


Perfect guidance!




Diamond instruments with guide pin – set 4384A

Even the most practiced eye and steady hand benefit from expert support. Here it is – in the shape of Komet diamond instruments with a 0.5 mm long guide pin to give you extra control when preparing crown margins. Obstructed view on the operation site is no longer a problem, for the controlled guidance of the instrument guarantees that the defined distance to the periodontium is kept.

These instruments are the ultimate choice to ensure a homogeneous cutting depth and to avoid an over-preparation. After only one use, you are going to wonder how you ever did without our diamond instruments with guide pin for the preparation of crowns. Contact your local Komet representative.



 HENRY SCHEIN®

Henry Schein Halas · 44 O'Dea Avenue · Waterloo NSW 2017 · Ph 1300 65 88 22 · Fax 1300 65 88 10 · info@henryschein.com.au · www.henryschein.com.au
 Henry Schein Shalfoon · 23 William Pickering Drive · Albany, North Shore · Auckland · Ph +64 9 414 0040 · Fax +64 9 414 0041 · info@henryschein.co.nz · www.henryscheinshalfoon.co.nz